



Apogee Prescription

To prescribe the Apogee, complete this form or use your own. Follow the same procedure as when ordering any portable oxygen prescription.

Fax the completed Rx to 636-778-1050.

- Patient Name: _____
- Patient Phone #: () -
- Patient Date of Birth: ____ / ____ / _____
- Applicable Diagnosis: _____
- Apogee Portable Oxygen System @ ____ flow setting (1 - 6) via nasal cannula
 - Or adjust flow setting to keep O₂ Saturation to ____ % via nasal cannula
- Dispense portable oxygen cylinder (circle):

M6 D E Other: _____
- Physician Name: _____
- Physician Signature: _____
- Date: ____ / ____ / _____
- NPI # : _____
- Office # : _____



Office: 636-778-1926
Fax: 636-778-1050



info@dynaris.com
www.dynaris.com



743 Spirit 40 Park Drive, Suite 108
Chesterfield, MO 63005



The Apogee System Duration by Cylinder Size and Setting

| Setting |  Apogee System with B (M6) Cylinder <i>4.9 lbs 164 Liters</i> |  Apogee System with D Cylinder <i>8.1 lbs 425 Liters</i> |  Apogee System with E Cylinder <i>10.7 lbs 680 Liters</i> |
|----------------|---|--|---|
| 1 | 13.5 hours | 35 hours | 56.5 hours |
| 2 | 7 hours | 17.5 hours | 28 hours |
| 3 | 4.5 hours | 11.5 hours | 18.5 hours |
| 4 | 3.5 hours | 8.5 hours | 14 hours |
| 5 | 2.5 hours | 7 hours | 11 hours |
| 6 | 2 hours | 5 hours | 9 hours |

Notes:

- Data calculated from the industry standard 20 breaths per minute.
- Times (in hours) are approximate.
- Weights (in lbs) include all necessary components of the Apogee System and cylinder.